**Instruction Manual
and Experiment Guide
for the PASCO scientific
Model SE-9638**

012-03471
11/87

INTRODUCTION

EQUIPMENT

OPERATION

Measuring e/m

Analysis of e/m Measurement

Deflections of Electrons

Two Simple Demonstrations

DIAGRAM OF CONNECTIONS

Page

1

1

2

2

3


4

4

5

e/m

APPARATUS

Copyright © November 1987

UNITED STATES \$5.00

TABLE OF CONTENTS

The PASCO Model 5B-5000 electron apparatus provides a simple method for measuring electron charge-to-mass ratio of the electron. The method is similar to that used by J.J. Thomson in 1897. A beam of electrons is accelerated through

vacuum. A pair of Helmholtz coils produces a uniform and measurable magnetic field at right angles to the electron beam. The deflection of the electron beam in a circular path is determined by the accelerating potential (V), the current in the coils (i), and the radius of the circular path (r). The charge-to-mass ratio (e/m) is easily calculated when these quantities are explained in the Operation section.

INTRODUCTION

EQUIPMENT

OPERATION

Measuring e/m

Analysis of e/m Measurement

Deflections of Electrons in an Electric Field

Two Simple Demonstrations

DIAGRAM OF CONNECTIONS

A magnetometer allows the electron beam to be deflected at any angle (from 0 to 180 degrees) with respect to the magnetic field from the Helmholtz coils. The magnetometer also allows you to examine the vector nature of the magnetic force by making charged particles. Other experiments are also possible with


Figure 1 The 5B-5000 electron apparatus

the magnetometer. The magnetometer also allows you to examine the effect of a magnetic field on the electron beam.

EQUIPMENT

The 5B-5000 electron apparatus (see Figure 1) is filled with helium gas at a pressure of 10^{-5} mm Hg, and contains an electron gun and deflection coils. The electron beam leaves a visible trail in the tube, because some of the electrons collide with the helium atoms, which are excited and then radiate visible light.

The electron gun is shown in Figure 2. The heater heats the cathode, which emits electrons. The electrons are accelerated by a potential applied between the cathode and the anode. The grid is held positively with respect to the cathode and repels electrons with respect to the anode. It helps to focus the electron beam.

CAUTION: The voltage on the heater of the electron gun should never exceed 6.3 volts. Higher voltages will burn out the heater and destroy the tube.

Figure 2 Electron gun

INTRODUCTION

The PASCO Model SE-9638 e/m apparatus provides a simple method for measuring e/m , the charge to mass ratio of the electron. The method is similar to that used by J.J. Thomson in 1897. A beam of electrons is accelerated through a known potential, so the velocity of the electrons is known. A pair of Helmholtz coils produces a uniform and measurable magnetic field at right angles to the electron beam. This magnetic field deflects the electron beam in a circular path. By measuring the accelerating potential (V), the current to the Helmholtz coils (I), and the radius of the circular path of the electron beam (r), e/m is easily calculated: $e/m = 2V/B^2r^2$. (The calculations are explained in the Operation section of this manual.)

The e/m apparatus also has deflection plates that can be used to demonstrate the effect of an electric field on the electron beam. This can be used as a confirmation of the negative charge of the electron and also to demonstrate how an oscilloscope works.

A unique feature of the e/m tube is that the socket rotates, allowing the electron beam to be oriented at any angle (from 0-90 degrees) with respect to the magnetic field from the Helmholtz coils. You can therefore rotate the tube and examine the vector nature of the magnetic forces on moving charged particles. Other experiments are also possible with


Figure 1 The e/m Apparatus

the e/m tube. For example, you can use a small permanent magnet instead of the Helmholtz coils to investigate the effect of a magnetic field on the electron beam.

EQUIPMENT

The e/m Tube — The e/m tube (see Figure 2) is filled with helium at a pressure of 10^{-2} mm Hg, and contains an electron gun and deflection plates. The electron beam leaves a visible trail in the tube, because some of the electrons collide with the helium atoms, which are excited and then radiate visible light.

The electron gun is shown in Figure 3. The heater heats the cathode, which emits electrons. The electrons are accelerated by a potential applied between the cathode and the anode. The grid is held positive with respect to the cathode and negative with respect to the anode. It helps to focus the electron beam.

CAUTION: The voltage to the heater of the electron gun should NEVER exceed 6.3 volts. Higher voltages will burn out the filament and destroy the e/m tube.


Figure 2 e/m Tube

The Helmholtz Coils—The geometry of Helmholtz coils—the radius of the coils is equal to their separation—provides a highly uniform magnetic field. The Helmholtz coils of the e/m apparatus have a radius and separation of 15 cm. Each coil has 130 turns. The magnetic field (B) produced by the coils is proportional to the current through the coils (I) times 7.80×10^{-4} tesla/ampere [B (tesla) = $(7.80 \times 10^{-4})I$].

The Controls—The control panel of the e/m apparatus is straightforward. All connections are labeled. The hook-ups and operation are explained in the next section.

Cloth Hood—The hood can be placed over the top of the e/m apparatus so the experiment can be performed in a lighted room.

Mirrored Scale—A mirrored scale is attached to the back of the rear Helmholtz coil. It is illuminated by lights that light automatically when the heater of the electron gun is powered. By lining the electron beam up with its image in the mirrored scale, you can measure the radius of the beam path without parallax error.

Additional Equipment Needed —

Power Supplies:

6-9 VDC @ 3 a (ripple < 1%) for Helmholtz coils


Figure 3 Electron Gun

(PASCO Model SF-9584 Low Voltage Power Supply)

6.3 VDC or VAC for filament

150-300 VDC accelerating potential (PASCO Model SF-9585 High Voltage Power Supply)

Meters:

Ammeter with 0-2 A range to measure current in Helmholtz coils (PASCO Model SE-9590 or SE-9545 Handheld Multimeter)

Voltmeter with 0-300 V range to measure accelerating potential (PASCO Model SE-9590 or SE-9545 Handheld Multimeter)

OPERATION

Measuring e/m

1. If you will be working in a lighted room, place the hood over the e/m apparatus.
2. Flip the toggle switch up to the e/m MEASURE position.
3. Turn the current adjust knob for the Helmholtz coils to the OFF position.
4. Connect your power supplies and meters to the front panel of the e/m apparatus, as shown in Figure 4.
5. Adjust the power supplies to the following levels:

CAUTION: The voltage to the heater of the electron gun should NEVER exceed 6.3 volts. High voltages will burn out the filament and destroy the e/m tube.

ELECTRON GUN

HEATER:	6.3 VAC or VDC
ELECTRODES:	150 to 300 VDC
HELMHOLTZ COILS:	6-9 VDC (ripple should be less than 1%)

6. Slowly turn the current adjust knob for the Helmholtz coils clockwise. Watch the ammeter and take care that the current does not exceed 2 A.
7. Wait several minutes for the cathode to heat up. When it does, you will see the electron beam emerge from the electron gun and it will be curved by the field from the Helmholtz coils. Check that the electron beam is parallel to the Helmholtz coils. If it is not, turn the tube until it is. Don't take it out of its socket. As you rotate the tube, the socket will turn.
8. Carefully read the current to the Helmholtz coils from your ammeter and the accelerating voltage from your voltmeter. Record the values below.

Current to Helmholtz coils = $I =$ _____

Accelerating voltage = $V =$ _____

9. Carefully measure the radius of the electron beam. Look through the tube at the electron beam. To avoid parallax errors, move your head to align the electron beam with the reflection of the beam that you can see


Figure 4 Connections for e/m Experiment

on the mirrored scale. Measure the radius of the beam as you see it on both sides of the scale, then average the results. Record your results below.

Electron Beam Radius = $r =$ _____

Analysis of e/m Measurement

The magnetic force (F_m) acting on a charged particle of charge q moving with velocity v in a magnetic field (B) is given by the equation $F_m = qv \times B$, (where F , v , and B are vectors and \times is a vector cross product). Since the electron beam in this experiment is perpendicular to the magnetic field, the equation can be written in scalar form as:

$$F_m = evB \quad (1)$$

where e is the charge of the electron.

Since the electrons are moving in a circle, they must be experiencing a centripetal force of magnitude

$$F_c = mv^2/r \quad (2)$$

where m is the mass of the electron, v is its velocity, and r is the radius of the circular motion. Since the only force acting on the electrons is that caused by the magnetic field, $F_m = F_c$, so equations 1 and 2 can be combined to give $evB = mv^2/r$ or

$$e/m = v/Br \quad (3)$$

Therefore, in order to determine e/m , it is only necessary to know the velocity of the electrons, the magnetic field produced by the Helmholtz coils, and the radius of the electron beam.

The electrons are accelerated through the accelerating potential V , gaining kinetic energy equal to their charge times the accelerating potential. Therefore $eV = 1/2 mv^2$, the velocity of the electrons is therefore:

$$v = (2eV/m)^{1/2} \quad (4)$$

The magnetic field produced near the axis of a pair of Helmholtz coils is given by the equation:

$$B = \frac{[N\mu_0]I}{(5/4)^{3/2}a} \quad (5)$$

A derivation for this formula can be found in most introductory texts on electricity and magnetism.

Equations 4 and 5 can be plugged into equation 3 to get a final formula for e/m :

$$e/m = v/Br = \frac{2V(5/4)^{3/2}a^2}{(N\mu_0 Ir)^2}$$

where:

V = the accelerating potential

a = the radius of the Helmholtz coils

N = the number of turns on each Helmholtz coil = 130

μ_0 = the permeability constant $= 4\pi \times 10^{-7}$

I = the current through the Helmholtz coils

r = the radius of the electrons beam path

Deflections of Electrons in an Electric Field

IMPORTANT: Do not leave the beam on for long periods of time in this mode. The beam will ultimately wear through the glass walls of the tube.

You can use the deflection plates to demonstrate how the electron beam is deflected in an electric field.

1. Setup the equipment as described above for measuring e/m except:
 - a. Flip the toggle switch to ELECTRICAL DEFLECT.
 - b. Do not supply current to the Helmholtz coils.
 - c. Connect a 0-50 VDC power supply between the banana plug connectors labeled DEFLECT PLATES (UPPER and LOWER)
2. Apply the 6.3 VDC or VAC to the HEATER and 150-300 VDC to the ELECTRODES of the ELECTRON GUN (the accelerating potential). Wait several minutes to warm up the cathode.

3. When the electron beam appears, slowly increase the voltage to the deflection plates from 0 V to approximately 50 VDC. Note the deflection of the electron beam. Note that the beam is bent towards the positively charged plate.

Two Simple Demonstrations

1. Instead of using the Helmholtz coils to bend the electron beam, you can use a permanent magnet to show the effect of a magnetic field on the electron beam. Just provide the following power to the e/m apparatus:

HEATER: 6.3 VAC or VDC
ELECTRON GUN ELECTRODES: 150-300 VDC

When the electron beam appears, use your permanent magnet to bend the beam.
2. The socket for the e/m tube is designed so that the tube can be rotated 90 degrees. The tube can therefore be oriented so it is at any angle, from 0-90 degrees, with respect to the magnetic field from the Helmholtz coils. By setting up the equipment as for measuring e/m, you can rotate the tube and study how the beam deflection it affected.


Diagram of Connections for the
SE-9638 e/m Apparatus

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PHYWE
excellence in science

A mass spectrometer

Mass spectrometry is an analytical technique, which accurately measures the mass of different molecules within a sample. It is normally used to identify samples and determine the purity of samples.

A mass spectrometer generates multiple ions from the sample under investigation, it then separates them according to their specific mass-to-charge ratio, and then records the relative abundance of each ion type. Results are displayed as spectra of the signal intensity of detected ions as a function of the mass-to-charge ratio.

Other information (1/2)

PHYWE
excellence in sciencePrior
knowledge

When electrons are bombarding into a gas sample, the gas break into charged molecules, or ionized. By subjecting them to an electric field and magnetic field, ions of the same mass-to-charge ratio will undergo the same amount of deflection.

Scientific
principle

Electrons are accelerated in an electric field and enter a magnetic field at right angles to the direction of motion. The specific charge of the electron is determined from the accelerating voltage, the magnetic field strength and the radius of the electron orbit.

Other information (2/2)

PHYWE
excellence in scienceLearning
objective

To learn about the specific charge of the electron, Lorentz force and electron ionization.

Tasks

Determination of the specific charge of the electron (e/m_0) from the path of an electron beam in crossed electric and magnetic fields of variable strength.

Safety instructions

PHYWE
excellence in science

For this experiment the general instructions for safe experimentation in science lessons apply.

At this experiment dangerous voltages are used. Under no circumstances wires and the plugs must be touched. Only the given at workstation high voltages shall be used. The heater voltages of the tubes to produce the electron beam shall not exceed the given voltages.

Be very careful when handling the structure.

Theory (1/3)

PHYWE
excellence in science

If an electron of mass m_0 and charge e is accelerated by a potential difference U , it attains the kinetic energy:

$$e \cdot U = \frac{1}{2} \cdot m_0 \cdot v^2$$

where v is the velocity of the electron.

In a magnetic field of strength \vec{B} , the Lorentz force acting on an electron with velocity \vec{v} is:

$$\vec{F} = e \cdot \vec{v} \times \vec{B}$$

Theory (2/3)

PHYWE
excellence in science

If the magnetic field is uniform, as it is in the Helmholtz arrangement, the electron therefore follows a spiral path along the magnetic lines of force, which becomes a circle of radius r if \vec{v} is perpendicular to \vec{B} .

Since the centrifugal force $m_0 \cdot v^2 / r$ thus produced is equal to the Lorentz force, we obtain

$$v = \frac{e}{m_0} \cdot B \cdot r$$

where B is the absolute magnitude of \vec{B} . From equation (1), it follows that

$$\frac{e}{m_0} = \frac{2U}{(Br)^2}$$

Theory (3/3)

PHYWE
excellence in science

To calculate the magnetic field B , the first and fourth Maxwell equations are used in the case where no time dependent electric fields exist. We obtain the magnetic field strength B_z on the z -axis of a circular current I for a symmetrical arrangements of two coils at a distance a from each other:

$$B_z = \mu_0 \cdot I \cdot R^2 + \{ (R^2 + (z - \frac{a}{2})^2)^{3/2} + (R^2 + (z + \frac{a}{2})^2)^{3/2} \}$$

with $\mu_0 = 1.257 \cdot 10^{-6} \frac{Vs}{Am}$ and R is radius of coil.

For the Helmholtz arrangement of two coils ($a = R$) with number of turns n in the centre between the coils one obtains

$$B = \left(\frac{4}{5}\right)^{3/2} \cdot \mu_0 \cdot n \frac{I}{R}$$


PHYWE
excellence in science

Setup and procedure

Setup (1/2)

PHYWE
excellence in science

Wiring diagram for Helmholtz coils

The two coils are turned towards each other in the Helmholtz arrangement. Since the current must be the same in both coils, connection in series is preferable to connection in parallel.


The maximum permissible continuous current of 5A should not be exceeded.

If the polarity of the magnetic field is correct, a curved luminous trajectory is visible in the darkened room.

Setup (2/2)

If the trace has the form of a helix this must be eliminated by rotating the narrow beam tube around its longitudinal axis.

For detailed description of the narrow beam tube, please refer to the operating instructions.

Wiring diagram for Narrow beam tube

Procedure

Experimental setup

By varying the magnetic field (current) and the velocity of the electrons (acceleration and focussing voltage), the radius of the orbit can be adjusted, such that it coincides with the radius defined by the luminous traces.

When the electron beam coincides with the luminous traces, only half of the circle is observable. The radius of the circle is then 2, 3, 4 or 5 cm.

Evaluation (1/2)

PHYWE
excellence in science

For the coils

$$R = 0.15 \text{ m}$$

$$n = 124$$

$$R = 0.2 \text{ m and } n = 154$$

current I and specific charge of the electron e/m_0 are determined for various voltages U and various radii r of the electron trajectories.

Compare the mean value of the measured specific charge of the electron and the literature value:

$$e/m = 1.759 \cdot 10^{11} \text{ As/kg}$$

		$r = 0.02 \text{ m}$		$r = 0.03 \text{ m}$		$r = 0.04 \text{ m}$		$r = 0.05 \text{ m}$	
$\frac{U}{V}$	I	$\frac{e/m_0}{10^{11} \frac{\text{As}}{\text{kg}}}$		$\frac{e/m_0}{10^{11} \frac{\text{As}}{\text{kg}}}$		$\frac{e/m_0}{10^{11} \frac{\text{As}}{\text{kg}}}$		$\frac{e/m_0}{10^{11} \frac{\text{As}}{\text{kg}}}$	
100	2.5	1.7		1.6	1.8	1.1	2.2	0.91	2.0
120	2.6	1.9		1.7	1.9	1.3	1.9	1.0	2.0
140	2.8	1.9		1.9	1.8	1.4	1.9	1.1	1.9
160	-	-		2.0	1.9	1.5	1.9	1.2	1.9

$$\overline{e/m_0} = (1.84 \pm 0.02) \cdot 10^{11} \text{ As/kg}$$

Current I and specific charge of the electron e/m_0

Evaluation (2/2)

PHYWE
excellence in science

How does Lorentz force acts on the electron?

☐ It causes the electron to move in a circular motion

☐ It influences the speed of the electron

☒ Check

What might be the most remarkable error in the measurements?

☐ Radius of the coils

☐ Electric current

☐ Radius of the circle

☒ Check